

# 27<sup>th</sup> Annual Scientific Congress

31 May - 2 June 2019

Level 2, Hong Kong Convention and Exhibition Centre



Hong Kong College of Cardiology

## Overcoming Challenges in LAAO: Looking Forward to a FLEXible Treatment Option

Dr Chui, Shing Fung  
Associate Consultant  
Queen Elizabeth Hospital  
Hong Kong



伊利沙伯醫院  
QUEEN ELIZABETH HOSPITAL


# Disclosure

Proctor      Boston Scientific Watchman Device

# Patient


- 78/F
- Hx of DM, HT, AF with CVA, previously on NOAC
- Recurrent massive GIB with extensive GI workup showed diverticulosis
- Referred for LAAO
- CHADS-VAS score 7
- HAS-BLED score 4

# Plan for a 27mm Watchman device



WL: 128 WW: 189

Cardiac  
Cardiac




Zoom: 278%  
Im: 1/59 Series: 2  
LittleEndianExplicit

1/3/2019, 10:07:16 AM  
Made In OsiriX

WL: 128 WW: 189

Cardiac  
Cardiac




Zoom: 278%  
Im: 1/55 Series: 3  
LittleEndianExplicit


1/3/2019, 10:16:05 AM  
Made In OsiriX

# Strategy



- ?Counter-clockwise rotation
- ?Deep implantation
- ?Kissing Watchman
- ?Use another device
  - Sandwich technique
- ?Continue oral anticoagulation



# Gentle “push” just before deployment




# PASS






# Incidence of AF increases with age





# Major Bleeding on OAC




No. at Risk

Apixaban	9088	8103	7564	5365	3048	1515
Warfarin	9052	7910	7335	5196	2956	1491


# LAA morphology




# Shallow LAA



# New Left Atrial Appendage Closure Device



**WATCHMAN FLX™**




**WATCHMAN™**

Designed to:

- Treat more patient anatomies, offer greater Flexibility
- Have improved maneuverability, ease-of-use
- Have enhanced stability, with new anchor design, and additional anchors
- Greater apposition to LAA wall, better seal and heal

# More Flexibility in Patient Treatment


Minimum LAA Depth =  $\frac{1}{2}$  Device Size



10-20%

**WATCHMAN FLX™**

Minimum LAA Depth  $\approx$  Ostium Dia.



**WATCHMAN™**

Less appendage depth needed for deployment


(i.e. 20mm device  $\approx$  10mm depth)

Polyethylene terephthalate fabric covering 90% of the device

# More Flexibility in Patient Treatment

Wider treatment range, greater device sizing overlap


## Sizing Charts



**WATCHMAN™**



**WATCHMAN FLX™**



**8-20%**  
Compression

**10-30%**  
Compression

Ostium Diameter (mm)

# Intuitive, safe and precise positioning

WATCHMAN FLX™ features a **soft, closed, atraumatic distal end** to enable advancement within the LAA when partially deployed.



The **distal fluoroscopic marker** allows high visibility of the device in relation to the LAA during device positioning and enhances atraumatic maneuverability.




# More Control and Ease of Use

Closed distal end allows greater maneuverability into the LAA

**WATCHMAN FLX™**



*FLX Ball*



Forms a bulb or ball when partially released, allowing safe antegrade negotiation of LAA depth

# More Control and Ease of Use

## WATCHMAN FLX™

'J' anchor




18 Anchors




It allows **Full or Partial Recapture**

## WATCHMAN™

'Straight' anchor



10 Anchors




Technically speaking,  
only allows **Partial Recapture**


# Better Seal and Heal

The 18 struts Nitinol frame conforms to even challenging LAA anatomies providing **optimal apposition of the long PET fabric to the LAA walls**, minimizing the risk of peri-device leaks.

The 18 anchors staggered in two rows ensure confident stability in the LAA.





# Threaded Insert Designs



**WATCHMAN™**

Exposed metal area: 6.7mm<sup>2</sup>




>70%  
reduction

**WATCHMAN FLX™**

Exposed metal area: 1.5mm<sup>2</sup>

Less  
Exposed  
Metal




Note: Values do not include the area of the threaded hole

# New Left Atrial Appendage Closure Device

**THREADED INSERT WITH  
REDUCED METAL  
EXPOSURE**

3,000 R&D HOURS



**OPTIMIZED FRAME  
SHAPE FOR SMOOTH  
DELIVERY**

7,500 R&D HOURS



**2 ROWS OF  
18 PRECISION  
ANCHORS**

24,000 R&D HOURS




**CONFORMABLE  
18-PAIR LASER  
CUT FRAME  
DESIGN**

10,000 R&D HOURS



**ATRAUMATIC  
CLOSED END**

10,000 R&D HOURS





# Deployment of WATCHMAN FLX™



FLX Ball Deployment Video



# Final Result of WATCHMAN FLX



# WATCHMAN'S Clinical Plan to Emcompasses

## All Patients with NVAF

Non-Valvular Atrial Fibrillation Patient Population





# PINNACLE FLX IDE STUDY

- Single arm non-randomized study design
  - DOAC only options for post-implant drug regimen
  - Non-inferiority to performance goal based on WATCHMAN 2.5
- Up to 490 enrollments (includes 90 roll-ins)
- Up to 45 US sites
- Primary Safety Endpoint
  - The occurrence of one of the following events between the time of implant and within 7 days following the procedure or by hospital discharge, whichever is later: all-cause death, ischemic stroke, systemic embolism, or device- or procedure- related events requiring open cardiac surgery or major endovascular intervention such as pseudoaneurysm repair, AV fistula repair, or other major endovascular repair.
- Primary Effectiveness Endpoint
  - The rate of effective LAA closure defined as any peri-device flow  $\leq$  5mm demonstrated by TEE at 12 months (US)
- Secondary Endpoints
  - The occurrence of ischemic stroke or systemic embolism at 24 months from the time of enrollment
- Follow-up at 45 days, 6, 12, 18 and 24 months



# Conclusion

- Incidence of AF increases with age
- Current generation LAA closure device still has limitations
- With modified, flexible designs, more patients with AF and high bleeding risks can be treated with the new Watchman FLX device

# Thank you

